


- Company
- Technology
- Applications
- Products
- Press & News
- Extranet
- Contact

Applications ▶▶

- Life Sciences
- Thin Films and Surface Treatment
  - Langmuir-Blodgett
  - Coatings
  - Crystallization
  - Thin film
  - Membranes
- Nanotubes and Nanowires
- Others

## Langmuir-Blodgett

### Gallery ▶▶



## Two-Component Phase Patterns in Model Pulmonary Surfactant Systems


**Collaboration: H. Nakahara & O. Shibata (Department of Biophysical Chemistry, Faculty of Pharmaceutical Sciences, Nagasaki International University, Japan)**

### Sarfus Benefits ▶▶

- ▶ Direct visualisation of LB layers at nanoscale
- ▶ Optical thickness measurement
- ▶ Large field of view

### Required Material ▶▶

- ▶ Standard surf
- ▶ Sarfus 3D
- ▶ Optical Microscope with DIC mode



### Comments ▶▶

Two-component monolayers (XHel 13-5 = 0.1) of dipalmitoylphosphatidylcholine (DPPC) and Hel 13-5 (or a mimicking peptide of SP-B protein in pulmonary surfactants) were transferred onto the Surf by the Langmuir-Blodgett (LB) technique [H. Nakahara et al., Langmuir 22 (2006) 1182-1192]. The LB film preparation was made at 35 mN/m on a 0.02 M Tris buffer and 0.13 M NaCl (pH 7.4) at 298.2 K. Considering the molecular thickness and orientation, the higher particle is made of DPPC and the surrounding lower region consists of Hel 13-5. The determined optical height difference of -0.7 and -1.1nm provides a good agreement with the previous data obtained by AFM measurements [H. Nakahara et al., Langmuir 22 (2006) 1182-1192].

ZAE La Pécardière  
Parc des Sittelles  
F - 72450 Montfort-le-Gesnois

Legal notice - Site map - Contact